



## Product Datasheet

### PNKD antibody (orb541038)

|                            |                                                        |
|----------------------------|--------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to PNKD                     |
| <b>Reactivity</b>          | Human, Mouse                                           |
| <b>Conjugation</b>         | Unconjugated                                           |
| <b>Tested Applications</b> | ICC, IF, WB                                            |
| <b>Immunogen</b>           | Recombinant fusion protein of human PNKD(NP_056303.3). |
| <b>Preservatives</b>       | PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.     |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.            |
| <b>Note</b>                | For research use only                                  |
| <b>Isotype</b>             | IgG                                                    |
| <b>Clonality</b>           | Polyclonal                                             |
| <b>MW</b>                  | 43kDa                                                  |
| <b>Uniprot ID</b>          | <b>Q8N490</b>                                          |
| <b>Dilution Range</b>      | WB: 1:500-1:2000, IF: 1:50-1:200                       |
| <b>Expiration Date</b>     | 12 months from date of receipt.                        |



Immunofluorescence analysis of U2OS cell...



Western blot analysis of extracts of var...