


Product Datasheet

PPP1R15A antibody (orb541096)

Description	Rabbit polyclonal antibody to PPP1R15A
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human PPP1R15A(NP_055145.3).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	110kDa
Uniprot ID	Q75807
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot
analysis of
extracts of
var...