

Product Datasheet

PREP antibody (orb541122)

Description	Rabbit polyclonal antibody to PREP
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human PREP(NP_002717.3).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	81kDa
Uniprot ID	P48147
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.