



Product Datasheet

PRKAA1 antibody (orb541126)


Description	Rabbit polyclonal antibody to PRKAA1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthetic peptide of human AMPK α 1(NP_006242.5).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	64kDa
Uniprot ID	Q13131
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human to...



Immunohistochemical staining of rat sple...



Western blot analysis of extracts of var...