

PRKCA antibody (orb541134)

Western blot analysis showing PRVc protein expression in various tissues. The lanes are labeled: HkA-1, HkA-2, HkA-3, HkA-4, HkA-5, HkA-6, HkA-7, HkA-8, HkA-9, HkA-10, HkA-11, HkA-12, HkA-13, HkA-14, HkA-15, HkA-16, HkA-17, HkA-18, HkA-19, HkA-20, HkA-21, HkA-22, HkA-23, HkA-24, HkA-25, HkA-26, HkA-27, HkA-28, HkA-29, HkA-30, HkA-31, HkA-32, HkA-33, HkA-34, HkA-35, HkA-36, HkA-37, HkA-38, HkA-39, HkA-40, HkA-41, HkA-42, HkA-43, HkA-44, HkA-45, HkA-46, HkA-47, HkA-48, HkA-49, HkA-50, HkA-51, HkA-52, HkA-53, HkA-54, HkA-55, HkA-56, HkA-57, HkA-58, HkA-59, HkA-60, HkA-61, HkA-62, HkA-63, HkA-64, HkA-65, HkA-66, HkA-67, HkA-68, HkA-69, HkA-70, HkA-71, HkA-72, HkA-73, HkA-74, HkA-75, HkA-76, HkA-77, HkA-78, HkA-79, HkA-80, HkA-81, HkA-82, HkA-83, HkA-84, HkA-85, HkA-86, HkA-87, HkA-88, HkA-89, HkA-90, HkA-91, HkA-92, HkA-93, HkA-94, HkA-95, HkA-96, HkA-97, HkA-98, HkA-99, HkA-100. The molecular weight markers are indicated on the left: 170kDa, 130kDa, 100kDa, 70kDa, 55kDa, 40kDa, 30kDa. A band labeled PRVc is visible in the 70kDa region.

Western blot
analysis of
extracts of
var...