



Product Datasheet

PSMA4 antibody (orb541164)


Description	Rabbit polyclonal antibody to PSMA4
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human PSMA4(NP_001096137.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	29kDa
Uniprot ID	P25789
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U2OS cell...



Immunohistochemical staining of human CO...



Western blot analysis of extracts of var...