



Product Datasheet

PSMB2 antibody (orb541168)


Description	Rabbit polyclonal antibody to PSMB2
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, IP, WB
Immunogen	Recombinant fusion protein of human PSMB2(NP_002785.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	23kDa
Uniprot ID	P49721
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of MCF-7 cells.



Immunohistochemical staining of human colon tissue.



Western blot analysis of extracts of various cell lines.