



Product Datasheet

PSMB7 antibody (orb541170)


Description	Rabbit polyclonal antibody to PSMB7
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human PSMB7(NP_002790.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	30kDa
Uniprot ID	Q99436
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of NIH-3T3 C...



Immunofluorescence analysis of C6 cells



Western blot analysis of extracts of var...