


Product Datasheet

PSPC1 antibody (orb541181)

Description	Rabbit polyclonal antibody to PSPC1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	Recombinant fusion protein of human PSPC1(NP_001035879.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	65kDa
Uniprot ID	Q8WXF1
Dilution Range	WB: 1:500-1:2000, IP: 1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...