


Product Datasheet

PTPRD antibody (orb541207)

Description	Rabbit polyclonal antibody to PTPRD
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human PTPRD(NP_002830.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	180-215kDa
Uniprot ID	P23468
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:100-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...