



Product Datasheet

RAB7A antibody (orb541248)


Description	Rabbit polyclonal antibody to RAB7A
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human RAB7(NP_004628.4).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	23KDa
Uniprot ID	P51149
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U-2OS cel...



Immunofluorescence analysis of NIH/3T3 C...



Immunofluorescence analysis of C6 cells

...