



Product Datasheet

RABIF antibody (orb541253)

Description	Rabbit polyclonal antibody to RABIF
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human RABIF(NP_002862.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	14kDa
Uniprot ID	P47224
Dilution Range	WB: 1:500-1:2000, IF 1:50 - 1:100
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of L929 cell...



Western blot analysis of extracts of var...