


Product Datasheet

RAD54B antibody (orb541261)

Description	Rabbit polyclonal antibody to RAD54B
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, IP, WB
Immunogen	Recombinant fusion protein of human RAD54B(NP_036547.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	103kDa
Uniprot ID	Q9Y620
Dilution Range	WB: 1:500-1:2000, IF: 1:50-1:200, IP: 1:50-1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...