


Product Datasheet

RAP2B antibody (orb541273)

Description	Rabbit polyclonal antibody to RAP2B
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human RAP2B(NP_002877.2).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	18kDa
Uniprot ID	P61225
Dilution Range	WB: 1:100-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...