


Product Datasheet

RBM15 antibody (orb541300)

Description	Rabbit polyclonal antibody to RBM15
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human RBM15(NP_001188474.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	120kDa
Uniprot ID	Q96T37
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot
analysis of
extracts of
var...