


Product Datasheet

RFPL2 antibody (orb541334)

Description	Rabbit polyclonal antibody to RFPL2
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human RFPL2(NP_006596.2).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	42kDa
Uniprot ID	O75678
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot
analysis of
extracts of
var...