


## Product Datasheet

### RNASE12 antibody (orb541366)

|                            |                                                       |
|----------------------------|-------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to RNASE12                 |
| <b>Reactivity</b>          | Human                                                 |
| <b>Conjugation</b>         | Unconjugated                                          |
| <b>Tested Applications</b> | WB                                                    |
| <b>Immunogen</b>           | A synthetic peptide of human RNASE12(NP_001019993.1). |
| <b>Preservatives</b>       | PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.      |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.           |
| <b>Note</b>                | For research use only                                 |
| <b>Isotype</b>             | IgG                                                   |
| <b>Clonality</b>           | Polyclonal                                            |
| <b>MW</b>                  | 17kDa                                                 |
| <b>Uniprot ID</b>          | <b>Q5GAN4</b>                                         |
| <b>Dilution Range</b>      | WB: 1:500-1:2000                                      |
| <b>Expiration Date</b>     | 12 months from date of receipt.                       |



Western blot analysis of extracts of var...