



Product Datasheet

ROBO1 antibody (orb541401)


Description	Rabbit polyclonal antibody to ROBO1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human ROBO1(NP_002932.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	270kDa
Uniprot ID	Q9Y6N7
Dilution Range	WB: 1:200-1:2000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U-2 OS cells...



Western blot analysis of extracts of various sources...



Western blot analysis of extracts of various sources...