

RPS26 antibody (orb541441)

Western blot analysis of RPS26 expression in various tissues. The blot shows bands for RPS26 (approximately 20 kDa) and a loading control (approximately 15 kDa). The lanes are labeled: U-531MG, UO2, Puji, DU145, Mouse brain, Mouse heart, Mouse kidney, and Mouse liver. Molecular weight markers are indicated on the left: 70kDa, 55kDa, 40kDa, 35kDa, 25kDa, 15kDa, and 10kDa. The RPS26 band is present in all lanes, with varying intensities.

Western blot
analysis of
extracts of
var...