



Product Datasheet

SCARB2 antibody (orb541502)

Description	Rabbit polyclonal antibody to SCARB2
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Recombinant fusion protein of human SR-B2/LIMPII(NP_005497.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	37kDa/54kDa
Uniprot ID	Q14108
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat heart...



Western blot analysis of extracts of var...