



## Product Datasheet

### SERPINE2 antibody (orb541566)

|                            |                                                            |
|----------------------------|------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to SERPINE2                     |
| <b>Reactivity</b>          | Mouse                                                      |
| <b>Conjugation</b>         | Unconjugated                                               |
| <b>Tested Applications</b> | WB                                                         |
| <b>Immunogen</b>           | Recombinant fusion protein of human SERPINE2(NP_006207.1). |
| <b>Preservatives</b>       | PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.         |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.                |
| <b>Note</b>                | For research use only                                      |
| <b>Isotype</b>             | IgG                                                        |
| <b>Clonality</b>           | Polyclonal                                                 |
| <b>MW</b>                  | 44kDa                                                      |
| <b>Uniprot ID</b>          | <b>P07093</b>                                              |
| <b>Dilution Range</b>      | WB: 1:500-1:2000, IHC-P: 1:50-1:200                        |
| <b>Expiration Date</b>     | 12 months from date of receipt.                            |



Immunohistochemical staining of mouse kidney tissue.



Western blot analysis of extracts of various mouse tissues.