



Product Datasheet

SERPINI1 antibody (orb541570)

Description	Rabbit polyclonal antibody to SERPINI1
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Recombinant fusion protein of human SERPINI1 (NP_005016.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	46kDa
Uniprot ID	Q99574
Entrez	5274
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:100, IF: 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat brain...



Western blot analysis of extracts of var...