


Product Datasheet

SETD2 antibody (orb541575)

Description	Rabbit polyclonal antibody to SETD2
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human SETD2(NP_054878.5).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	260kDa
Uniprot ID	Q9BYW2
Dilution Range	WB: 1:200-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...