



Product Datasheet

SFRP1 antibody (orb541582)

Description	Rabbit polyclonal antibody to SFRP1
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human SFRP1(NP_003003.3).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	35kDa
Uniprot ID	Q8N474
Dilution Range	WB: 1:100-1:500, IHC-P: 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human st...



Western blot analysis of extracts of var...