


## Product Datasheet

### SLC17A6 antibody (orb541638)

|                            |                                                    |
|----------------------------|----------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to SLC17A6              |
| <b>Reactivity</b>          | Mouse, Rat                                         |
| <b>Conjugation</b>         | Unconjugated                                       |
| <b>Tested Applications</b> | WB                                                 |
| <b>Immunogen</b>           | A synthetic peptide of human SLC17A6(NP_065079.1). |
| <b>Preservatives</b>       | PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.   |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.        |
| <b>Note</b>                | For research use only                              |
| <b>Isotype</b>             | IgG                                                |
| <b>Clonality</b>           | Polyclonal                                         |
| <b>MW</b>                  | 64kDa                                              |
| <b>Uniprot ID</b>          | <b>Q9P2U8</b>                                      |
| <b>Dilution Range</b>      | WB: 1:500-1:2000                                   |
| <b>Expiration Date</b>     | 12 months from date of receipt.                    |



Western blot analysis of extracts of var...