



Product Datasheet

SLC1A1 antibody (orb541641)


Description	Rabbit polyclonal antibody to SLC1A1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human SLC1A1(NP_004161.4).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	70KDa
Uniprot ID	P43005
Dilution Range	WB: 1:500-1:2000, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of mouse bra...



Immunofluorescence analysis of rat brain...



Western blot analysis of extracts of var...