



Product Datasheet

SLC2A2 antibody (orb541668)

Description	Rabbit polyclonal antibody to SLC2A2
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human GLUT2/GLUT2/SLC2A2(NP_000331.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	38-45KDa
Uniprot ID	P11168
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human li...



Western blot analysis of extracts of var...