



Product Datasheet

SLC2A3 antibody (orb541669)

Description	Rabbit polyclonal antibody to SLC2A3
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human GLUT3/SLC2A3(NP_008862.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	48-60KDa
Uniprot ID	P11169
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse he...



Western blot analysis of extracts of var...