

Product Datasheet

SLC01C1 antibody (orb541718)

Description	Rabbit polyclonal antibody to SLC01C1
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human SLC01C1(Q9NYB5).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	78kDa
Uniprot ID	Q9NYB5
Dilution Range	WB: 1:200-1:2000
Expiration Date	12 months from date of receipt.