



Product Datasheet

SPTAN1 antibody (orb541788)


Description	Rabbit polyclonal antibody to SPTAN1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human Alpha Fodrin(NP_001123910.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	285kDa
Uniprot ID	Q13813
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U-2 OS cells...



Immunofluorescence analysis of U-2 OS cells...



Western blot analysis of extracts of various tissues...