


Product Datasheet

SRSF7 antibody (orb541802)

Description	Rabbit polyclonal antibody to SRSF7
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human SRSF7(NP_001026854).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	35KDa
Uniprot ID	Q16629
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.



Western blot
analysis of
extracts of
var...