



Product Datasheet

SSTR2 antibody (orb541815)


Description	Rabbit polyclonal antibody to SSTR2
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	A synthetic peptide of human Somatostatin Receptor 2 (SSTR2)(NP_001041.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	56kDa
Uniprot ID	P30874
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U-2 OS cells



Immunofluorescence analysis of L929 cells



Immunofluorescence analysis of C6 cells