



Product Datasheet

ST3GAL4 antibody (orb541823)

Description	Rabbit polyclonal antibody to ST3GAL4
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human ST3GAL4(NP_006269.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	34-38kDa
Uniprot ID	Q11206
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat kidney tissue.



Western blot analysis of extracts of various tissues.