



Product Datasheet

STIP1 antibody (orb541842)

Description	Rabbit polyclonal antibody to STIP1
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, IP, WB
Immunogen	Recombinant fusion protein of human STIP1(NP_006810.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	69kDa
Uniprot ID	P31948
Dilution Range	WB: 1:500-1:2000, IF: 1:50-1:200, IP: 1:20-1:50
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...



Western blot analysis of extracts of var...