


Product Datasheet

SUCLG1 antibody (orb541860)

Description	Rabbit polyclonal antibody to SUCLG1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human SUCLG1(NP_003840.2).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	26-36kDa
Uniprot ID	P53597
Dilution Range	WB: 1:200-1:2000
Expiration Date	12 months from date of receipt.



Western blot
analysis of
extracts of
var...