



Product Datasheet

SUCLG2 antibody (orb541861)

Description	Rabbit polyclonal antibody to SUCLG2
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human SUCLG2(NP_003839.2).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	47KDa
Uniprot ID	Q96I99
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.




Immunofluorescence
analysis of U-2 OS
ce...



Immunofluorescence analysis of C6 cells

...



Western blot
analysis of extracts
of var...