



Product Datasheet

SUV39H1 antibody (orb541881)

Description	Rabbit polyclonal antibody to SUV39H1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human SUV39H1(NP_003164.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	48KDa
Uniprot ID	O43463
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF: 1:50-1:200, ChIP: 1:20-1:100
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of C6 cells



Western blot analysis of extracts of var...