


Product Datasheet

TACR2 antibody (orb541898)

Description	Rabbit polyclonal antibody to TACR2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human TACR2 (NP_001048.2).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	44kDa
Uniprot ID	P21452
Entrez	6865
Dilution Range	WB: 1:200-1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...