

Product Datasheet

TADA3 antibody (orb541899)

Description	Rabbit polyclonal antibody to TADA3
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human TADA3(NP_597814.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	55kDa
Uniprot ID	075528
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.

