


Product Datasheet

TAS1R3 antibody (orb541918)

Description	Rabbit polyclonal antibody to TAS1R3
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human TAS1R3(NP_689414.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	93kDa
Uniprot ID	Q7RTX0
Dilution Range	WB: 1:1000-1:4000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...