

TBRG4 antibody (orb541934)

Western blot analysis showing TBR levels in various tissues. The lanes are labeled: OVCAR3, H1-29, U-251MG, 293T, Mouse spleen, Mouse brain, Mouse heart, Mouse kidney, and Rat heart. Molecular weight markers are indicated on the left: 170kDa, 130kDa, 100kDa, 70kDa, 55kDa, 40kDa, and 25kDa. TBR bands are visible in the 70kDa and 55kDa regions across all lanes.

Western blot
analysis of
extracts of
var...