

Product Datasheet

TBXA2R antibody (orb541937)

Description	Rabbit polyclonal antibody to TBXA2R
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human TBXA2R(NP_001051.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	45kDa
Uniprot ID	P21731
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.