


Product Datasheet

TDO2 antibody (orb541955)

Description	Rabbit polyclonal antibody to TDO2
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human TDO2(NP_005642.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	48kDa
Uniprot ID	P48775
Dilution Range	WB: 1:200-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...