


Product Datasheet

TEKT5 antibody (orb541961)

Description	Rabbit polyclonal antibody to TEKT5
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human TEKT5(NP_653275.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	56kDa
Uniprot ID	Q96M29
Dilution Range	WB: 1:200-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...