

Product Datasheet

TFR2 antibody (orb541977)

Description	Rabbit polyclonal antibody to TFR2
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human TFR2(NP_003218.2).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	100KDa
Uniprot ID	Q9UP52
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.