



Product Datasheet

TGF beta 1 antibody (orb541979)


Description	Rabbit Polyclonal antibody to TGF beta 1
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF
Immunogen	Recombinant fusion protein of human TGF beta 1(NP_000651.3).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	44kDa
Uniprot ID	P01137
Dilution Range	WB: 1:500-1:2000, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U2OS cell...



Immunofluorescence analysis of NIH/3T3 C...



Immunofluorescence analysis of C6 cells

...