



Product Datasheet

TGF beta Receptor 2 antibody (orb541982)


Description	Rabbit Polyclonal antibody to TGF beta Receptor 2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Recombinant fusion protein of human TGF beta Receptor II (TGF beta Receptor II (TGFR2)) (NP_003233.4).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	90kDa
Uniprot ID	P37173
Entrez	7048
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human ma...



Immunohistochemical staining of human br...



Western blot analysis of extracts of var...