



Product Datasheet

TGF beta Receptor 2 antibody (orb541983)


Description	Rabbit Polyclonal antibody to TGF beta Receptor 2
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human TGF beta Receptor II (TGFB2)(NP_003233.4).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	90KDa
Uniprot ID	P37173
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U2OS cell...



Immunohistochemical staining of human ma...



Immunohistochemical staining of human br...