


Product Datasheet

TIMM50 antibody (orb542009)

Description	Rabbit polyclonal antibody to TIMM50
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human TIMM50(NP_001001563.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	39kDa
Uniprot ID	Q3ZCQ8
Dilution Range	WB: 1:1000-1:4000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...