


## Product Datasheet

### TNFRSF13C antibody (orb542063)

|                            |                                                             |
|----------------------------|-------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to TNFRSF13C                     |
| <b>Reactivity</b>          | Human, Mouse                                                |
| <b>Conjugation</b>         | Unconjugated                                                |
| <b>Tested Applications</b> | WB                                                          |
| <b>Immunogen</b>           | Recombinant fusion protein of human TNFRSF13C(NP_443177.1). |
| <b>Preservatives</b>       | PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.            |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.                 |
| <b>Note</b>                | For research use only                                       |
| <b>Isotype</b>             | IgG                                                         |
| <b>Clonality</b>           | Polyclonal                                                  |
| <b>MW</b>                  | 19kDa                                                       |
| <b>Uniprot ID</b>          | <a href="https://www.uniprot.org/entry/Q96RJ3">Q96RJ3</a>   |
| <b>Dilution Range</b>      | WB: 1:500-1:2000                                            |
| <b>Expiration Date</b>     | 12 months from date of receipt.                             |



Western blot analysis of extracts of var...